

## Introducing

# Topi-CLICK® MICRO™



Revolutionary Way to Accurately & Easily Micro-Dose  
Topical Medications for Humans & Pets



Targeted Topical  
Micro-Dosing

Easy Screw Cap



Choice of three  
applicators




Standard Luer  
Lock Connection



REFILL Indicator



Innovative  
TRIAD Shape  
Allows Precision  
Application  
Control



Three Visual  
CLICK Points  
Each turn/CLICK of  
Base dispenses 0.05 mL  
of topical

1 CLICK = 0.05 mL  
2 CLICKs = 0.10 mL  
3 CLICKs = 0.15 mL

## The Topi-CLICK® Micro™ Difference

### Accurate & Consistent Dose™

- Topi-CLICK® dispensers proven most accurate among others tested by \*ARL
- Eliminates user-variability
- Consistent automated dosing throughout treatment

### Automated-Metered Dosing

- Audible, Tactile and Visual CLICKs indicate measure of dose at each point of the unique triad-shaped base (One CLICK = 0.05 mL)

### Targeted Topical Therapy

- Innovative triad-shaped dispenser allows for precision handling during application
- Applicator reduces the risk of hand transference to others
- Offers applicator selections for customized dosing

Select one of 3  
applicators



Soft  
Angled  
Applicator™



Rounded  
Applicator™



Pin Point  
Applicator™

### Viscosity Assurance™

- Accuracy NOT affected by varying formulations
- Accuracy NOT affected by varying bases

Dispense/CLICK: 0.05 mL

Dispenser Volume: Up to 9.75 mL

Total Delivery Volume: Up to 9.50 mL

Approx. Residual: 0.15 mL

Approx. Priming: 0.10 or 2 CLICKs

Made of FDA grade materials | Phthalate free | blocks UV and is Recyclable ♻️.

## Topi-CLICK® MICRO™




### Topical Micro-dosing Application

- ✓ Automated-metered topical dosing provides an accurate and consistent dose every time
- ✓ Audible, Tactile and Visual CLICKs indicate measure of dose (no need to manually mark dose)
- ✓ Innovative TRIAD shape allows for precision handling
- ✓ Built-in REFILL Badge helps pet owners reorder medication on-time



**STEP 1:**  
Unscrew the cap of the Topi-CLICK® Micro™ and set aside for recapping after use.



**STEP 2:**  
Holding the Micro™ in an upright position, turn/CLICK the **Base** the number of CLICKs or mLs specified on the prescription label. Use the Audible, Tactile, and Visual CLICK points of the triad shaped **Base** as dosing indicators. Be sure to align **Base** with the **Triad Tube** to ensure dose is fully dispensed.



**STEP 3:**  
Apply topical medication with built-in applicator as directed by veterinarian or pharmacist. **Continue rubbing until preparation has completely saturated the prescribed areas through your pet's fur.** **IMPORTANT:** The built-in applicator lessens the risk of transference to unintended areas. No need to use your fingertips to apply medication.



**STEP 4:**  
Screw on cap after each dose to prevent evaporation. Be sure cap is secured completely.

Innovative TRIAD Shape Allows Precision Application Control

Each turn/Click of the base dispenses 0.05 mL of medication accurately.

### REFILL Badge™



REFILL Line

Plunger

**The Refill Badge™ reminder system helps pet owners reorder in a timely manner.**


When the **TOP** black line of the **Plunger** reaches the embossed line below the word "REFILL" approximately 32 CLICKs are left. Approximately 12 CLICKs remain at the line above the word "REFILL".

*Note: Follow instructions on prescription label before use. Topi-CLICK® Micro™ is not re-fillable.*

Made of FDA grade materials | Phthalate free | blocks UV and is Recyclable ♻️.

## Topi-CLICK® MICRO™

**Micro Quick-Fill Technology™**  
Makes Filling Fast & Simple



**Note: Graduated Lines** are for filling only and not for patients. **Top** black line of the plunger indicates **Fill Volume**. The embossed **Refill Badge** helps patients reorder medication in a timely manner.



Select one of 3 applicators



**Note: Soft Angled Applicator™** placement requires tip to be aligned where lower slant of tip is on the same side as **Refill Badge of the Micro™ tube**.

**Step 1: Filling**  
**Syringe Tip-to-Tip Filling-**  
Screw Luer syringe to the Luer lock tip of the Micro™. **Requires no extra Luer Lock adaptor.**

**EMP Jar Quick-fill Adaptor™ Filling-** Screw on Topi-CLICK® Micro™ EMP Jar Quick-fill Adaptor™. Fill topical medications directly from mixing jar with this added Micro™ accessory.

**Step 2: Applicator Tip Attachment -**  
Select Applicator Tip and press on by pushing straight down (important) against a hard surface.



**Step 3: Priming the Micro** Advance the actuator to engage plunger by manually turning/CLICKing base. The Actuator Bar should seat completely to the plunger before priming. **Actuator can also be advanced by using a mini cordless screwdriver to speed up the process (200 RPM Max). A Micro Driver Socket Bit™ can be purchased for faster actuator advancement with a cordless mini screwdriver.** (See Advancing Actuator Instructions.)




**Step 4: Cap Placement**  
Screw on cap securely to protect tip and prevent evaporation.



**Step 5: Label Application**  
Apply above the base opposite the "REFILL Badge" indicator.

**Step 3a - Priming**



Base

**Step 3a: Priming**  
Prime the Micro™ by turning/CLICKing the base until contents are dispensed, then wipe off. Be sure to align Micro™ base with the triad tube.  
**Approximate Priming: 0.10 mL or 2 CLICKS.**

**Dispense/CLICK: 0.05 mL**  
**Dispenser Volume: Up to 9.75 mL**  
**Total Delivery Volume: Up to 9.50 mL**

**Approx. Residual: 0.15 mL**  
**Approx. Priming: 0.10 or 2 CLICKS**



Micro Quick-Fill Technology<sup>™</sup>  
Makes Filling Fast & Simple

## Instructions to Advance the Micro<sup>™</sup> Actuator with a Mini Cordless Screwdriver



### Advancing the Actuator with a Mini Cordless Screwdriver

**Actuator** can be advanced by using a mini cordless screwdriver without a reversible trigger (200 RPM Max). Accidentally turning the base backwards could cause the Micro to slip in your hand.

The BLACK & DECKER<sup>™</sup> 3.6 Volt Pink Rechargeable Electric Screwdriver is recommended.

A **Micro Driver Socket Bit<sup>™</sup>** can also be purchased for faster actuator advancement with a cordless screwdriver. (See info below.)

### Topi-CLICK<sup>®</sup> Micro Driver Socket Bit<sup>™</sup> - SKU#: TC0910

The triad-shaped Micro driver is designed to perfectly conform to the Micro Base. It should be used with a 1/4 inch hex bit power screw driver. We recommend power drivers with a speed of 200 RPM or less to reduce the likelihood of damage to the Micro.

**Step 1:** Insert the **Micro Driver** into the screwdriver socket until it is properly seated. Depending on the manufacturer and model, there is usually a light click that lets you know when it is properly seated in the tool.


**Step 2:** Rotate the **Driver** in a clockwise direction to advance the **Micro Actuator** up the tube.



### Notes:

*We recommend stopping the Actuator Bar 1 mL short of the target graduation (currently only for priming 6 mL or smaller prescriptions.)*

**Step 3:** Run up the **Actuator** using the **Actuator Bar** as a visual indicator. The **Actuator Bar** should seat completely to the plunger.







## DoseLogix Laboratory Test Results

Product performed per design specifications. Device is both highly accurate and precise, demonstrating a tight variance and an average dose being 98.7% of the rated dose. Product performance variance was reduced for two CLICKs versus one CLICK.

Qualifies inspection and is passed.

### Accurate Micro-dosing Designed for Patients and Pets

#### Topi-CLICK® Micro™ Lifetime Dosing Accuracy

(1 CLICK or 0.05 mL)

(2 CLICKS or 0.10 mL)

Average Dose	0.049 mL	0.099 mL
Standard Deviation	0.003 mL	0.004 mL
Coefficient of Variation	5.7%	3.8%
Average Dose to Rated Volume	98.7%	98.7%

Dispense/CLICK: 0.05 mL  
Dispenser Volume: Up to 9.75 mL  
Total Delivery Volume: 9.50 mL

Approx. Residual: 0.15 mL  
Approx. Priming: 0.10 mL or 2 CLICKS

\*Volumes are approximate. As required for compounded prescription pharmaceuticals dispensed in volume dispensers; verify the weight to volume conversions thoroughly.

Made of FDA grade materials | Phthalate free | blocks UV and is Recyclable ♻️