

NEW!

Industry advancement in topical micro-dosing for humans and pets

Topi-CLICK[®]
MICRO[™]



Choose Topi-CLICK[®] Micro[™] for The Best User Experience[™]



✓ Accurate metered micro-dosing for patients

✓ Choice of 3 applicator styles for customized dosing



Rounded Applicator[™]



Soft Angled Applicator[™]



Pin Point Applicator[™]

✓ Micro Quick-fill Technology[™] for compounders

Topi-CLICK[®] Micro[™] Patent Pending

Micro Quick-fill Technology[™]
Save time and money in the pharmacy

Fill directly from a syringe with Micro's[™] Luer Lock Tip



Fill directly from a mixing jar with the Micro EMP Jar Quick-fill Adaptor[™]



ONE is easier than a **HANDFUL!**




Request **FREE** Topi-CLICK[®] Micro[™] Samples at
DoseLogix.com | 877.870.8448 | info@DoseLogix.com

DOSELOGIX
Copyright © 2018 DoseLogix, LLC All Rights Reserved.

Introducing

Topi-CLICK® MICRO™



Targeted Topical
Micro-Dosing

Revolutionary way to accurately & easily micro-dose topical medications for humans & pets

Easy Screw Cap



Choice of three applicators



Standard Luer Lock Connection



REFILL Indicator



The Topi-CLICK® Micro™ Difference

Accurate & Consistent Dose™

- Topi-CLICK® Micro™ is the most accurate among those in a patient simulation test.
- Eliminates user-variability
- Consistent automated dosing throughout treatment

Automated Metered Dosing

- Audible, Tactile and Visual CLICKs indicate measure of dose at each point of the unique triad-shaped base (One CLICK = 0.05 mL)

Targeted Topical Therapy

- Innovative triad-shaped dispenser allows for precision handling during application
- Applicator reduces the risk of hand transference to others
- Offers applicator selections for customized dosing

Select one of 3 applicators



Rounded Applicator™



Soft Angled Applicator™




Pin Point Applicator™

Viscosity Assurance™

- Accuracy NOT affected by varying formulations
- Accuracy NOT affected by varying bases

Innovative TRIAD Shape Allows Precision Application Control



Each turn/CLICK of Base dispenses 0.05 mL of topical

1 CLICK = 0.05 mL
2 CLICKs = 0.10 mL
3 CLICKs = 0.15 mL

Dispense/CLICK: 0.05 mL
Dispenser Volume: Up to 9.75 mL
Total Delivery Volume: Up to 9.45 mL

Approx. Residual: 0.15 mL
Approx. Priming: 0.15 or 3 CLICKs

Made of FDA grade materials | Phthalate free | blocks UV and is Recyclable ♻️.


TOPI-CLICK[®] MICRO[™] User Instructions



- ✓ Automated, metered topical dosing provides an accurate and consistent dose every time
- ✓ Audible, Tactile and Visual CLICKs indicate measure of dose (no need to manually mark dose)
- ✓ Innovative triad-shape allows for precision handling
- ✓ Built-in REFILL Badge helps patients reorder medication on-time



STEP 1:
Unscrew the cap of the Topi-CLICK[®] Micro[™] and set aside for recapping after use.



STEP 2:
Holding the Micro[™] in an upright position, turn/CLICK the **Base** the number of CLICKs or mLs specified on the prescription label. Use the Audible, Tactile and Visual CLICK points of the triad-shaped **Base** as dosing indicators. Be sure to align **Base** with the **Triad Tube** to ensure dose is fully dispensed.



STEP 3:
Apply topical medication with built-in applicator as directed by physician or pharmacist. **Continue rubbing until preparation is completely absorbed into the skin.**
IMPORTANT: The built-in applicator lessens the risk of transference to unintended areas. No need to use your fingertips to apply medication.



STEP 4:
Screw on cap after each dose to prevent evaporation. Be sure cap is secured completely.



Each turn/Click of the base dispenses 0.05 mL of medication accurately.



The Refill Badge[™] reminder system helps patients reorder in a timely manner.

When the **TOP** black line of the **Plunger** reaches the embossed line below the word "REFILL" approximately 32 CLICKs are left. Approximately 12 CLICKs remain at the line above the word "REFILL".

Note: Graduated Marks are for pharmacy filling purposes only. Follow instructions on prescription label before use. Topi-CLICK[®] Micro[™] is not re-fillable.

Made of FDA grade materials | Phthalate free | blocks UV and is Recyclable ♻️.

TOPI-CLICK® MICRO™

Micro Quick-Fill Technology™
Makes Filling Fast & Simple



Note: Graduated Lines are for filling only and not for patients. **Top** black line of the plunger indicates **Fill Volume**. The embossed **Refill Badge** helps patients reorder medication in a timely manner.



Select one of 3 applicators



Note: Soft Angled Applicator™ placement requires tip to be aligned where lower slant of tip is on the same side as **Refill Badge of the Micro™ tube**.

Step 1: Filling
Syringe Tip-to-Tip Filling-
Screw Luer syringe to the Luer lock tip of the Micro™. **Requires no extra Luer Lock adaptor.**

EMP Jar Quick-fill Adaptor™ Filling- Screw on Topi-CLICK® Micro™ EMP Jar Quick-fill Adaptor™. Fill topical medications directly from mixing jar with this added Micro™ accessory.

Step 2: Applicator Tip Attachment -
Select Applicator Tip and press on by pushing straight down (important) against a hard surface.



Step 3: Priming the Micro Advance the actuator to engage plunger by manually turning/CLICKing base. The Actuator Bar should seat completely to the plunger before priming. **Actuator can also be advanced by using a mini cordless screwdriver to speed up the process (200 RPM Max). A Micro Driver Socket Bit™ can be purchased for faster actuator advancement with a cordless mini screwdriver.** (See Advancing Actuator Instructions.)




Step 4: Cap Placement
Screw on cap securely to protect tip and prevent evaporation.



Step 5: Label Application
Apply above the base opposite the "REFILL Badge" indicator.

Step 3a - Priming



Base

Step 3a: Priming
Prime the Micro™ by turning/CLICKing the base until contents are dispensed, then wipe off. Be sure to align Micro™ base with the triad tube.
Approximate Priming: 0.15 mL or 3 CLICKs.

Dispense/CLICK: 0.05 mL
Dispenser Volume: Up to 9.75 mL
Total Delivery Volume: Up to 9.45 mL

Approx. Residual: 0.15 mL
Approx. Priming: 0.15 mL or 3 CLICKs



TOPI-CLICK[®]
MICRO[™] Micro Quick-Fill Technology[™]
Makes Filling Fast & Simple

Instructions to Advance the Micro[™] Actuator with a Mini Cordless Screwdriver



Advancing the Actuator with a Mini Cordless Screwdriver

Actuator can be advanced by using a mini cordless screwdriver without a reversible trigger (200 RPM Max). Accidentally turning the base backwards could cause the Micro to slip in your hand.

The BLACK & DECKER[™] 3.6 Volt Pink Rechargeable Electric Screwdriver is recommended.

A **Micro Driver Socket Bit**[™] can also be purchased for faster actuator advancement with a cordless screwdriver. (See info below.)

Topi-CLICK[®] Micro Driver Socket Bit[™] - SKU#: TC0910

The triad-shaped Micro driver is designed to perfectly conform to the Micro Base. It should be used with a 1/4 inch hex bit power screw driver. We recommend power drivers with a speed of 200 RPM or less to reduce the likelihood of damage to the Micro.

Step 1: Insert the **Micro Driver** into the screwdriver socket until it is properly seated. Depending on the manufacturer and model, there is usually a light click that lets you know when it is properly seated in the tool.


Step 2: Rotate the **Driver** in a clockwise direction to advance the **Micro Actuator** up the tube.



Notes:

We recommend stopping the Actuator Bar 1 mL short of the target graduation (currently only for priming 6 mL or smaller prescriptions.)

Step 3: Run up the **Actuator** using the **Actuator Bar** as a visual indicator. The **Actuator Bar** should seat completely to the plunger.





DoseLogix Laboratory Test Results

Product performed per design specifications. Device is both highly accurate and precise, demonstrating a tight variance and an average dose being 98.7% of the rated dose. Product performance variance was reduced for two CLICKs versus one CLICK.

Qualifies inspection and is passed.

Accurate Micro-dosing Designed for Patients and Pets

Topi-CLICK® Micro™ Lifetime Dosing Accuracy

(1 CLICK or 0.05 mL)

(2 CLICKS or 0.10 mL)

Average Dose	0.049 mL	0.099 mL
Standard Deviation	0.003 mL	0.004 mL
Coefficient of Variation	5.7%	3.8%
Average Dose to Rated Volume	98.7%	98.7%

Dispense/CLICK: 0.05 mL
Dispenser Volume: Up to 9.75 mL
Total Delivery Volume: Up to 9.45 mL

Approx. Residual: 0.15 mL
Approx. Priming: 0.15 mL or 3 CLICKS

*Volumes are approximate. As required for compounded prescription pharmaceuticals dispensed in volume dispensers; verify the weight to volume conversions thoroughly.

Made of FDA grade materials | Phthalate free | blocks UV and is Recyclable