



Topi-CLICK
by DOSELOGIX

**Industry Leader in
Metered-dosing Applicators**

Proud to Support
Compounding Pharmacists
as You Serve Your
Patients & Communities

**ACCURATE & EASY
Dosing Solutions for
Humans & Pets**

Innovative Vaginal Dosing

NEW
Rose
Quartz
Color!

#1
in Accuracy!

0.25 mL
per
CLICK

*American Journal
of Internal Medicine*
Perl® ARL Study Evaluation

Topi-CLICK
Perl

Micro-dosing for Pets

0.05 mL
per
CLICK

Topi-CLICK
MICRO

#1
in Accuracy!

0.25 mL
per
CLICK

Topi-CLICK
35

**Topi-CLICK® DoseCheck™
for Healthcare Facility Use**

Features
graduations
to record
dose.

0.25 mL
per
CLICK

Now Available
with 3Port
Applicator Pad!

Topi-CLICK 35
DOSE ✓

Micro-dosing for Humans

0.05 mL
per
CLICK

Topi-CLICK
MICRO

Oral Dosing for Pets

0.25 mL
per
CLICK

Topi-CLICK
35

TO ORDER, CONTACT ONE OF OUR VALUED DISTRIBUTORS



Topi-CLICK Micro® patent pending.

DOSELOGIX
Copyright © 2020, DoseLogix. All rights reserved.

Request the *American Journal of Internal Medicine* Perl® ARL Study Evaluation,
the Topi-CLICK® ARL Accuracy Study and FREE Samples for Patient Trials at -
www.DoseLogix.com | 877.870.8448 | info@DoseLogix.com

Proven #1 in Accuracy!¹



Topi-CLICK® Perl VAGINAL DOSING



NEW! Breakthrough Revealed in
Feminine Dosing Study

Pharmacists - Dispense medications in the Topi-CLICK Perl®
- based on the *American Journal of Internal Medicine*
Perl® ARL Study Evaluation Results.

"The Topi-CLICK Perl® system performed with a **high level of accuracy and precision** in delivering doses of feminine cream, outperforming the other tested devices in all measures. "

"Evaluation of Topical Dispenser Systems for Vaginal Delivery of Feminine Medications."
American Journal of Internal Medicine, 2019; p.99.

Dosing Accuracy of each Vaginal Dosing Method



36% less medication wasted than
traditional feminine dosing methods.
"Single-Blinded Feminine Dosing Study." ARL BioPharma, Inc. 2019; p.19.



0.25 mL
per
CLICK

- Accurate
- Easy-to-Use
- Elegant Design

NEW COLOR



Topi-CLICK
Perl® in
Rose Quartz

Multiple color options
make it easy for patients
to distinguish between
their medications.


DOSELOGIX

Copyright © 2020, DoseLogix. All rights reserved.

DoseLogix.com | 877.870.8448 | info@DoseLogix.com

¹Request samples and get the complete
Topi-CLICK Perl® ARL Accuracy Study at
DoseLogix.com/samples/



Topi-CLICK® 35 DOSE ™
Topical Dosing Dispenser

Healthcare providers and caregivers can easily verify a dose with the
NEW graduated Topi-CLICK® 35 DoseCheck™ system.




*Provides an independent
check system to verify dose.*



Features graduated
marks at each 0.5 mL.

BOTTOM of Plunger
indicates the amount of
medication remaining.

Each turn/CLICK of base
dispenses 0.25 mL of topical
1 CLICK = 0.25 mL
2 CLICKs = 0.50 mL

The Topi-CLICK® 35 DOSE ™ topical dosing dispenser provides the same advanced technology as the Topi-CLICK® 35 and offers EXTRA assurance with graduated marks so the dose starting and stopping points can be easily verified and recorded. *Topi-CLICK® 35 DoseCheck™ initial start point should begin when the bottom of the plunger is at the 30 mL graduated mark.*

Consistent & Accurate Dose™

- Topi-CLICK® dispensers proven most accurate among others tested by *ARL
- Eliminates user-variability
- Consistent automated dosing throughout treatment
- Allows for easy recording of numerical volumes dispensed

Built-in Applicator Pad

- Reduces the risk of cross-contamination to others

CLICK-Assist™ Technology

- Clear Audible CLICKs
- Built-in Visual CLICK points
- Strong Tactile Sensations at each CLICK point

Viscosity Assurance™

- Accuracy NOT affected by varying formulations
- Accuracy NOT affected by varying bases

Made of FDA grade materials | BPA and Phthalate free | UV block and is Recyclable 

Request

*ARL PEER REVIEWED EVIDENCE at
http://doselogix.com/accuracy_study.



Copyright © 2020, DoseLogix. All rights reserved.

Request FREE Samples for Patient Trials at -
www.DoseLogix.com | 877.870.8448 | info@DoseLogix.com



Provides an independent check system to verify dose.

Healthcare Provider Use Instructions

STEP 1:

Healthcare provider should begin by reading dose volume while holding the dispenser at eye level.

Record volume by reading the graduated line at the **bottom of the plunger** and then record by using the healthcare facility's preferred method of dose recording. **(Figure a.)**

Holding dispenser above or below eye level may cause errors in reading. Topi-CLICK[®] 35 DoseCheck™ initial start point should begin when the bottom of the plunger is at the 30 mL graduated mark.



STEP 2:

Lift off the cap of the Topi-CLICK[®] dispenser by placing your thumb under the tab on the protective outer cap and pressing upward.




STEP 3:

Watch the visual indicators when rotating base. Be sure to leave the base in the registered CLICK position by aligning the base column with the tube node after each use.



STEP 4:

Turn the base to dispense the number of CLICKs or mLs specified in the instructions on the prescription label.



STEP 5:

Apply topical medication as directed by the physician or pharmacist. Rub the area of the skin using the **built-in applicator pad**. Continue rubbing until topical is completely absorbed into the skin. **IMPORTANT:** The built-in applicator lessens the risk of cross-contamination and ensures the full dose is applied to intended area.




STEP 6:

Replace cap after each dose to prevent evaporation. Be sure cap snaps on tightly.



STEP 7:

Read volume stopping point by holding dispenser at eye level and recording as indicated by your healthcare facility. (Special Note: See **Figure b.**)



Follow instructions on prescription label before using the Topi-CLICK[®].

Micro-dosing for Humans



Micro-dosing for Pets



**Topi-CLICK Micro[®]
is NOW available in
11 Vibrant Colors!**



**Discover the cost
savings of using
Topi-CLICK Micro[®]
over 1 mL syringes.**



✓ Accurate metered micro-dosing
for patients

✓ Choice of 3 applicator styles for
customized dosing



Rounded
Applicator™



Soft Angled
Applicator™



Pin Point
Applicator™

✓ Micro Quick-fill Technology™ for
compounders

TO ORDER, CONTACT ONE OF OUR VALUED DISTRIBUTORS



B&B Pharmaceuticals Inc.



Topi-CLICK Micro[®] patent pending.

Topi-CLICK
MICRO[®]

Choose ONE topical micro-dosing applicator

"The time and easiness of loading Micro[®] was a definite selling point and the feedback has all been positive. People like that they don't have to deal with numerous 1 mL syringes and leakage."

Nessa B. Denham, CPhT, Access Pharmacy



\$12.00

**Bottom Line Immediate Cost Reduction When Using
Topi-CLICK Micro[®] Instead of 1 mL Syringes
for a 90 Day Prescription**

\$10.00

Micro[®] offers
a \$4.26
Cost Reduction

\$8.00

\$6.00

\$4.00

\$2.00

\$0.00



**\$5.33
total
pharmacy
cost**



Topi-CLICK
MICRO[®]

**\$9.59
total
pharmacy
cost**



(9) 1 mL Syringes

Topi-CLICK Micro[®] is the most accurate among those tested in a patient simulation study.

Topi-CLICK Micro[®] patent pending.

DOSELOGIX


Copyright © 2020, DoseLogix. All rights reserved.

Request FREE Samples for Patient Trials at -
www.DoseLogix.com | 877.870.8448 | info@DoseLogix.com



Value Cost Comparison of Topi-CLICK Micro® to 1 mL Syringes for a 90 Day Prescription

Topical micro-dosing for humans and pets



Topi-CLICK Micro®		Times Qty. Needed	Cost	1 mL Syringes		Times Qty. Needed	Cost
Total Dispenser + Accessories							
Dispenser Cost	\$4.83	1	\$4.83	\$.25	9	\$2.25	
Labels (\$.08/each)	\$.08	2 (1 for Micro and 1 for patient bag)	\$.16	Labels (\$.08/each)	10 (9 for syringes and 1 for patient bag)	\$.80	
Bags (\$.08/each)	\$.08	1	\$.08	Bags (\$.08/each)	1	\$.08	
Tip Caps (\$.08/each)	n/a	n/a	n/a	Tip Caps (\$.08/each)	9	\$.72	
Finger Cots (\$.05/each)	n/a	n/a	n/a	Finger Cots (\$.05/each)	90	\$4.50	
Total Micro™ Dispenser + Accessories			\$5.07	Total 1 mL Syringes + Accessories			\$8.35
Prep Labor (\$15/hr. or .0042/sec)							
Filling Time	30 sec	1	\$.13	15 sec	9	\$.57	
Labeling Time	15 sec	2	\$.13	15 sec	10	\$.63	
Bagging Time	1 sec	1	\$.004	1 sec	9	\$.04	
Total Labor Micro®			\$.26	Total Labor 1 mL Syringes			\$1.24
Micro® Dispenser + Accessories + Prep Labor			\$5.33	1 mL Syringes Dispenser + Accessories +Prep Labor			\$9.59

Topi-CLICK Micro® is the most accurate among those tested in a patient simulation study.

Topi-CLICK Micro® patent pending.



Copyright © 2020, DoseLogix. All rights reserved.

Request FREE Samples for Patient Trials at -
www.DoseLogix.com | 877.870.8448 | info@DoseLogix.com