

NEW Topi-CLICK® UD™ 45



Topi-CLICK®
by Cordica Medical

An Optimized Experience and Improved Bottom Line

The **Topi-CLICK® UD™ 45 Metered-Dose Topical Applicator** is designed to enhance pharmacy efficiency and patient ease of use.

Convenient New Size!

Patients get a 3-month supply when dosing at 0.50 mL (two CLICKs) per day.

Patient, Physician and Pharmacist **Benefits**

- 01/ Supports **consistent, repeatable** dosing
- 02/ Enables **clean, controlled application** of topical dosage to intended area
- 03/ Integrated applicator **reduces the risk of contamination** compared to traditional containers
- 04/ Allows **verification of dispensed volume** through visible graduations and high-contrast indicator line
- 05/ **Audible, visual, and tactile** feedback during dispensing for improved dosing confidence
- 06/ Enables **in-device mixing**, reducing product transfer and handling during preparation
- 07/ Improves pharmacy **workflow efficiency** with integrated fast prime feature



NEW!

Topi-CLICK® UD™ 45

Metered-Dose Topical Applicator

This single-use device repeatably dispenses a metered volume of a topical cream or gel for application to the skin. It also functions as a container for compounding creams and gels with electronic mortar and pestle (EMP) systems. The Topi-CLICK® UD™ 45 applicator can be used for a variety of topically applied or transdermally absorbed medicament.

Product Specifications

Patented design, registered as an **FDA Class I** medical device

Dose size of 0.25 mL* per quarter-turn CLICK of base

Medicament contacting materials:
polyolefins free of phthalates and bisphenol A (BPA)

UV-inhibiting materials used for the medicament container

45 mL

Max usable volume

1.5 mL

Residual volume

**Independent accuracy study data available upon request*


Innovative Design Details



Protective cover securely locks into place for convenient storage and transport.



Applicator cap convex, five-port design facilitates even dispensing and hands-free application.



Threaded tube accepts mixing lids of most EMP systems, allowing cream and medicament to be mixed within the applicator.

Innovative **mechanical drive system** enables fast priming and precise dispensing.

Easy-to-read **graduated marks** for pharmacists, compounders and patients.

The **base** advances the plunger to dispense dosage, alerting the user by sound, sight and feel.

WHAT'S INCLUDED?

Each 25/pack includes a reusable prime tool. Mixing lids and additional paddles are available for purchase.